



FRENOVO

Dengue IgG/IgM/NS1 Combo Rapid Test Kit

Instructions For Use

PRODUCT NAME
FRENOVO Dengue IgG/IgM/NS1 Combo Rapid Test Kit
PACKAGESPECIFICATION
20 tests/kit
INTENDED USE

FRENOVO Dengue IgG/IgM/NS1 Combo Rapid Test is an in vitro qualitative Immuno chromatographic assay for the rapid detection of Dengue NS1 antigens and anti-dengue IgG/IgM antibodies in human blood, serum and plasma specimens simultaneously. The test results are intended to aid in the diagnosis of dengue infection.

SUMMARY AND PRINCIPLES OF THE PROCEDURE

Dengue fever is one of the most important mosquito-borne diseases in the world in the terms of morbidity, mortality. Dengue fever virus (serotypes 1 – 4) belongs to the group flavivirus, and is transmitted in nature by day-biting Aceder mosquitos. The most important mosquito vector is highly domesticated and urban species, Aedes aegypti. Primary Dengue infection, also known as Dengue Fever, is the most common type of dengue illness. It is associated with mild to high fever, headache, muscle pain and skin rash. Secondary infection is known as Dengue Hemorrhagic Fever (DHF) or Dengue Shock Syndrome, and often results in high fever and in many cases, with hemorrhagic events and circulatory failure. The fatality rate in patients with Dengue Shock Syndrome can be as high as 44%. Dengue presents typically as a fever of sudden onset with headache, retrobulbar pain, pain in the back and limbs (break-bone fever), lymphadenopathy and maculopapular rash. Patients diagnosed with dengue in endemic areas generally have secondary infection, whereas patients in non-endemic areas are usually diagnosed with primary infection. Specific antibody responses to Dengue virus enable serodiagnosis and differentiation between primary and secondary dengue infections.

Dengue NS1 (nonstructural protein I) is a highly conserved glycoprotein. NS1 antigen was found circulating in samples of infected patient from the first day up to 9 days after the onset of the fever. After the anti-NS1 antibody elevate in human body, the detectable NS1 antigens decline quickly. While detectable NS1 antigens declining, the elevated antibodies can be detected for longer time. Usually IgM does not become detectable until 3 to 10 days after the onset of illness in cases of primary dengue infection and until 2 to 3 days after onset of illness in secondary infections. In primary infections, IgG appear the 14th day and may persist for many years. Secondary infections generate an anamnestic IgG antibody response that is characterized by a rapid rise in IgG antibodies detectable at 4-5 days after the onset of the illness. The detection of both NS1 antigen and anti-Dengue antibodies provides the tool for the diagnosis of dengue infection from the early infection to the stage after the onset of the illness. It can enhance the accuracy of the diagnosis of dengue infection.

Dengue NS1 Test is a sandwich solid phase immunochromatographic assay. When sample is added to sample pad, it moves through the conjugate pad and mobilizes gold antiNS1 conjugate that is coated on the conjugate pad. The mixture moves along the membrane by capillary action and reacts with antiNS1 antibody that is coated on the test region. If NS1 is present, the result is the formation of a colored band in the test line region. If there is no NS1 in the sample, the area will remain colorless. The sample continues to move to the control area and forms a pink to purple color, indicating the test is working and the result is valid.

Dengue IgG/IgM Test utilizes the principle of Immunochromatography. Mouse anti human IgM and human IgG antibodies are immobilized on the nitrocellulose membrane respectively, as two individual test lines (IgM line and IgG line) in the test window of the test device. The IgG line in the test window is closer to the sample well and followed by IgM line. As the test sample flows through the membrane within the test device, the colored–Dengue specific recombinant antigen-colloidal gold conjugate complexes with specific antibodies (IgM and/or IgG) of Dengue virus, if present in the sample. This complex moves further on the membrane to the test region where it is captured by the antihuman IgM and/or human IgG antibodies coated on the membrane leading to formation of a colored band, which indicates a positive test results. Absence of this colored band in the test window indicates a negative test result. A built-in control line will always appear in the test window when the test has performed properly, regardless of the presence or absence of anti-Dengue virus antibodies in the specimen.

MATERIALS PROVIDED

1. Test cassettes: 20 pieces test cassettes individually pouched.
2. Wash Buffer Solution: 2 pieces dropper bottle with 2.5 ml solution.
3. Droppers: 20 pieces droppers of 25 ul for Dengue NS1 Test, 20 pieces droppers of 10 ul for Dengue IgG/IgM Antibody Test.
4. Package insert: 1 piece attached.

MATERIALS REQUIRED BUT NOT PROVIDED

- Timer or stopwatch.
- Specimen collection containers
- Disposable gloves and/or protective clothing
- Centrifuge(for plasma only)
- Micro-pipette

WARNINGS

1. Read the package insert completely before using the product. The instructions must be followed carefully as not doing so may result in inaccurate results.
2. The kits for diagnostic use only.
3. Perform test at room temperature.

PRECAUTIONS

1. FRENOVO Dengue IgG/IgM/NS1 Combo Rapid Test Kit is for professional use only.
2. The package insert instructions must be followed to ensure optimum test performance.
3. The kit is intended for in vitro diagnostic use.
4. As with all screening assays, any results should be considered presumptive until confirmatory assays have been performed according to local practice or WHO guidelines.

Safety Precautions

1. Standard precautions for handling infectious agents should be observed when using this kit.
2. Wear protective clothing such as lab coat, safety glasses and disposable gloves when handling specimens and assay reagents.

3. Wash hands thoroughly after use.
4. In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.

Bio safety Precautions

Appropriate bio safety practices should be used when handling specimens and reagents. These precautions include, but are not limited to the following:

1. Do not smoke ,eat, drink, apply cosmetics or handle contact lenses in areas in which specimens are handled.
2. Dispose of all specimens, used devices and tubes as though they are capable of transmitting infection. The preferred methods of disposal are by autoclave at 121°C for a minimum of 60 minutes or by incineration. Disposable materials may be incinerated. Liquid waste may be mixed with appropriate chemical disinfectants. A solution of 10% bleach is recommended. Allow 60 minutes for effective decontamination. NOTE: Do not autoclave solutions containing bleach.
3. When disposing of solution, avoid contact with acid to prevent liberation of a toxic gas.
4. All spills should be wiped thoroughly using a suitable disinfectant such as a sodium hypochlorite solution.
5. Use a separate dropper and device for each specimen tested.

Handling Precautions

1. Do not use if the kit safety seal is absent, damaged or broken.
2. Do not use any device if the pouches have been perforated.
3. Each device is for single use only.
4. Do not mix wash buffer solution/test cassettes from different kit lots.
5. Do not use the kit past the expiration date (this date is printed on the kit box).
6. Adequate lighting is required to read the test results.
7. The result should be read immediately after the end of the 10 minutes incubation time following the addition of specimen and wash buffer solution. Do not read results beyond 15 minutes.

STORAGE INSTRUCTIONS

1. The kit should be stored between 2-30°C and the shelf life is 24 months.
2. The kit components are stable until the expiration date printed on the outer label, when stored as directed. The kit expiry date is determined by whichever of the components has the shortest expiry date. The kit expiry date is not impacted once the wash buffer solution has been opened. Do not use kit components beyond overall kit expiry date.
3. If stored refrigerated, ensure that the pouched device is brought to room temperature before opening.
4. Do not freeze the kit.

SAMPLE COLLECTION AND PREPARATION

1. Applicable samples: Whole Blood/Serum/Plasma.
2. Separate serum or plasma from whole blood as soon as possible to avoid hemolysis. Use only clear nonhemolysis specimens.
3. Testing should be performed immediately after the specimens have been collected as soon as possible . Do not leave the specimens at room temperature for prolonged periods.
4. Serum and plasma specimens may be stored at 2–8 °C for up to 7 days, for long term storage,serum/plasma specimens should be kept below -20 °C. Whole blood collected by venipuncture should be stored at 2–8 °C if the test is to be run within 2 days of collection.
5. Do not freeze whole blood specimens.
6. Bring specimens to room temperature prior to testing. Frozen specimens must be completely thawed and mixed well prior to testing. Specimens should not be frozen and thawed repeatedly.
7. If specimens are to be shipped, they should be packed in compliance with local regulations covering the transportation of etiological agents.
8. EDTA K2, Heparin sodium, Citrate sodium and Potassium Oxalate can be used as the anticoagulant for collecting the specimen.

QUALITY CONTROL

An internal procedural control is included in the test. a colored line appearing in the control line region (C) is an internal valid procedural control, it confirming adequate membrane wicking. Control standards are not supplied with this kit; however, it is recommended that positive and negative controls be tested as a good laboratory practice to confirm the test procedure and to verify proper test performance.

TEST PROCEDURE

Allow the test cassette, specimen, buffer solution to equilibrate to room temperature (15-30°C) prior to testing.

1. Remove the test cassette from the sealed pouch and use it within one hour. Place the test cassette on a clean and level surface.
2. Dengue IgG/IgM Test

For Serum or Plasma Specimens


To use a dropper: Hold the dropper vertically, draw the specimen and transfer the specimen to the sample well of the test cassette (one drop/approximately 10ul), then squeeze the wash buffer solution bottle, add three drops wash buffer solution inside (approximately 90ul) to the sample well and start the timer. Avoid trapping air bubbles in the sample well.

To use a micro-pipette: Pipette and dispense 10ul of specimen to the sample well of the test cassette, then squeeze the wash buffer solution bottle, add three drops wash buffer solution inside (approximately 90ul) to the sample well and start the timer.

For Whole Blood Specimens

To use a dropper: Hold the dropper vertically, draw the specimen and transfer the specimen to the sample well of the test cassette (two drops/approximately 20ul), then squeeze the wash buffer solution bottle, add three drops wash buffer solution inside (approximately 90ul) to the sample well and start the timer. Avoid trapping air bubbles in the sample well.

To use a micro-pipette: Pipette and dispense 20ul of specimen to the sample well of the test cassette, then squeeze the wash buffer solution bottle, add three drops wash buffer solution inside (approximately 90ul) to the sample well and start the timer.



3. Dengue NS1 Test

For Serum or Plasma Specimens

To use a dropper: Hold the dropper vertically, draw the specimen and transfer the specimen to the sample well of the test cassette (two drops/approximately 50ul), then squeeze the wash buffer solution bottle, add one drop wash buffer solution inside (approximately 30ul) to the sample well and start the timer. Avoid trapping air bubbles in the sample well.



To use a micro-pipette: Pipette and dispense 50ul of specimen to the sample well of the test cassette, then squeeze the wash buffer solution bottle, add one drop wash buffer solution inside (approximately 30ul) to the sample well and start the timer.

For Whole Blood Specimens

To use a dropper: Hold the dropper vertically, draw the specimen and transfer the specimen to the sample well of the test cassette (two drops/approximately 50ul), then squeeze the wash buffer solution bottle, add one drop wash buffer solution inside (approximately 30ul) to the sample well and start the timer. Avoid trapping air bubbles in the sample well.

To use a micro-pipette: Pipette and dispense 50ul of specimen to the sample well of the test cassette, then squeeze the wash buffer solution bottle, add one drop wash buffer solution inside (approximately 30ul) to the sample well and start the timer.

4. Wait for the colored line(s) to appear. The test result should be read at 10 minutes. Do not interpret the result after 15 minutes.



INTERPRETATION OF RESULTS

Negative: One colored line should be in the control line region (C), no line detection lines appear, indicating that Dengue IgG/IgM Antibody or Dengue NS1 antigen have not been detected and the results are negative.

IgG POSITIVE: Two lines appear on Dengue IgG/IgM test strip. One colored line should be in the control line region (C), and a colored line appears in IgG test line region. The result is positive for Dengue virus specific-IgG and is probably indicative of secondary Dengue infection.

IgM POSITIVE: Two lines appear on Dengue IgG/IgM test strip. One colored line should be in the control line region (C), and a colored line appears in IgM test line region. The result is positive for Dengue virus specific-IgM antibodies and is indicative of primary Dengue infection.

IgG and IgM POSITIVE: Three lines appear. One colored line should be in the control line region (C), and two colored lines should appear in IgG test line region and IgM test line region. The color intensities of the lines do not have to match, the result is positive for IgG & IgM antibodies and is indicative of secondary Dengue infection.

Dengue NS1 Positive: Two lines appear on Dengue NS1 test strip, both of the quality control line C and the detection line T, the Dengue NS1 antigen have been detected and the result is positive.

NOTE: The intensity of the color in the IgG and/or IgM test line region(s) will vary depending on the concentration of Dengue antibodies in the specimen. Therefore, any shade of color in the IgG and/or IgM test line region(s) should be considered positive.

INVALID: Control line fails to appear. Insufficient buffer volume or incorrect procedural techniques are the most likely reasons for control line failure. Review the procedure and repeat the procedure with a new test cassette. If the problem persists, discontinue using the test kit immediately and contact your local distributor.

LIMITATIONS

- The test is for qualitative detection of Dengue NS1 antigen and/or anti-Dengue antibodies in whole blood, serum or plasma sample and does not indicate the quantity of the analytes.
- The test is for in vitro diagnostic use only.
- A negative test result cannot exclude a recent infection.
- It is common that viruses belong to flavivirus genus, such as dengue virus, Japanese encephalitis, tick-borne encephalitis, yellow fever virus and West Nile virus, have the serological cross reaction on antibody tests.
- As in case of all diagnostic tests, a definitive clinical diagnosis should not be based on the result of a single test but should rather be made after all the clinical findings have been evaluated.

PERFORMANCE CHARACTERISTICS

Sensitivity and Specificity

Clinical study was performed to compare the results obtained by The kit and PCR. The results indicated that The kit has a high sensitivity and specificity as summarized below:

Dengue IgM Study		PCR		Total Results
Dengue IgG/IgM test	Results	IgM Positive	IgM Negative	
	IgM Positive	249	9	
	IgM Negative	12	771	
Total Results		261	780	1041

Dengue IgM Study Summary Results:
Clinical sensitivity =95.40% (95%CI* 92.11% ~ 97.60%) Clinical specificity = 98.85% (95%CI* 97.82% ~ 99.47%)
Accuracy=97.98% (95%CI* 96.93% ~ 98.75%)

Dengue IgG Study		PCR		Total Results
Dengue IgG/IgM test	Results	IgG Positive	IgG Negative	
	IgG Positive	128	12	
	IgG Negative	12	889	
Total Results		140	901	1041

Dengue IgG Study Summary Results:
Clinical sensitivity =91.43% (95%CI* 85.51% ~ 97.60%) Clinical specificity = 98.85% (95%CI* 97.82% ~ 99.47%)
Accuracy=97.98% (95%CI* 96.93% ~ 98.75%)

Clinical study		PCR		Total Results
Dengue NS1 Test	Results	Positive	Negative	
	Positive	58	7	
	Negative	4	93	
Total Results		62	100	162

Accuracy Results:
Clinical sensitivity=93.55 % (95%CI* 84.30% to 98.21%) Clinical specificity=93.00 % (95%CI* 86.11% to 97.14%)
Accuracy=93.21% (95%CI* 88.18% to 96.56%)

Interference Substances

The following potential interfering substances have been tested using The kit and no interference was observed:

Dengue IgG/IgM Antibodies Interference Substance	Tested Concentration
Acetaminophen	136 µg/mL
Albumin	47 mg/mL
Bilirubin	7.40 mg/mL
Cholesterol	260 mg/mL
Creatinine	6.90 mg/mL
Digoxin	3.00 ng/mL
Ethanol	2.10 mg/mL
Glucose	3.90 mg/mL
Immunoglobulin A	2.20 mg/mL
Immunoglobulin G	10.00 mg/mL
Immunoglobulin M	1.10 mg/mL
Total Protein	70.40 mg/mL
Triglycerides	1.90 mg/mL
Urea nitrogen	700 µg/mL
Uric acid	92 µg/mL

Dengue NS1 antigen Interference Substance

Flavivirus	Strain	Concentration(TCID50)
St. Louis Encephalitis Virus	Parton	1*10 ^{5.82} U/mL
West Nile Virus	NY2001-6263	1*10 ^{5.66} U/mL
West Nile Virus 2	1986	1*10 ^{6.2} U/mL
West Nile Virus Lineage 2	B-956 Uganda	1*10 ^{5.18} U/mL
Yellow Fever Virus	17D	1*10 ^{5.02} U/mL
Zika Virus	N/A	1*10 ^{4.55} U/mL

Cross Reaction

Cross-reactivity of The kit was evaluated using serum samples containing antibodies to other pathogens., which has no effect on the negative and positive test results, and there is no cross-reaction.

Dengue IgM potential cross-reactant	Dengue IgG potential cross-reactant
Influenza A virus (H1N1, H3N2)	Influenza A virus (H1N1, H3N2)
Influenza B virus (Yamagata IgM,Victoria IgM)	Influenza B virus (Yamagata IgG, Victoria IgG)
Endemic human coronavirus (OC43, 229E)	Endemic human coronavirus (OC43, 229E)
CMV IgM	CMV IgG
Rubella IgM	Rubella IgG
Toxo IgM	Toxo IgG
HSV IgM	HSV IgG
Coxsackie virus group B IgM	Coxsackie virus group B IgG
Epstein-Barr virus IgM	Epstein-Barr virus IgG
Enterovirus 71 IgM	Enterovirus 71 IgG
Coxsackie virus type A16 IgM	Coxsackie virus type A16 IgG
Varicella zoster virus IgM	Varicella zoster virus IgG
Mumps Virus IgM	Mumps Virus IgG
Respiratory syncytial virus IgM	Respiratory syncytial virus IgG
Adenovirus IgM	Adenovirus IgG
Chlamydia pneumoniae IgM	Chlamydia pneumoniae IgG
Mycoplasma pneumoniae IgM	Mycoplasma pneumoniae IgG
Measles virus IgM	Measles virus IgG

Dengue NS1 potential cross-reactant	hepatitis C
measles	hepatitis B
rubella	leptospirosis
influenza A	hemorrhagic fever
epidemic cerebrospinal meningitis	sepsis
typhoid	tsutsugamushi
AIDS	Encephalitis B

Minimum detection limit

Dengue NS1 sub-type	minimum detection limit
Type I	8.3x10 ⁻⁷ TCID50/ml
Type II	8.3x10 ⁻⁷ TCID50/ml
Type III	1.5x10 ⁻⁷ TCID50/ml
Type IV	3.0x10 ⁻⁷ TCID50/ml

INDEX OF SYBOML

IVD	In vitro diagnostic medical device	ⓧ	single-use, Please don't reuse it
⌚	Use-by date	📖	Consult instructions for use
⚠️	Cautions	🏭	Manufacturer
🌡️ 30°C	Temperature limit	LOT	Batch code
📅	Date of manufacture	☔	Keep Dry
💡	Avoid overexposure to the sun	🚫	Don't use the product when the package is damaged
CE	CE mark	🦠	Biological risks

Hangzhou Frenovo Biotech Co., Ltd.
Address: Room 401, Building 26, Room 401&301, Building 27, No. 488-1, Donghu North Road, Donghu Street, LinPing District, Hangzhou City, Zhejiang Province, China.
Tel: 86-0571-89170657
Email: business@frenovo.com

CMC Medical Devices & Drugs S.L
Address: C/Horacio Lengo Nº 18, CP 29006, Málaga, Spain
Email: info@cmcmedicaldevices.com

INSTRUCTION APPROVAL AND REVISION DATE

Revision Date: 2023.01.23
Approval Date: 2023.01.23
Date of Issue: 2023.01.23