



**FRENOVO**

## Human Chorionic Gonadotropin (hCG) Rapid Test Kit Instructions For Use

### PRODUCT NAME

Human Chorionic Gonadotropin (hCG) Rapid Test Kit

### PACKAGE SPECIFICATION

20 tests/kit

### INTENDED USE

FRENOVO Human Chorionic Gonadotropin (hCG) Rapid Test is a lateral flow, qualitative immunoassay. It is intended for qualitative detection of hCG, with a minimum detection concentration of 25 mIU/mL in human urine samples, aim to assist of clinical diagnosis of pregnancy.

### SUMMARY AND PRINCIPLES OF THE PROCEDURE

Human chorionic gonadotropin (hCG) is a glycoprotein hormone produced by the developing placenta shortly after fertilization. In normal pregnancy, hCG can be detected in urine as early as 7 to 10 days after conception. 1,2,3,4 hCG levels continue to rise very rapidly, frequently exceeding 100 mIU/mL by the first missed menstrual period, 2,3,4 and peaking in the 100,000-200,000 mIU/mL range about 10-12 weeks into pregnancy. The appearance of hCG in urine soon after conception, and its subsequent rapid rise in concentration during early gestational growth, make it an excellent marker for the early detection of pregnancy.

FRENOVO Human Chorionic Gonadotropin (hCG) Rapid Test is a rapid test that qualitatively detects the presence of hCG at the sensitivity of 25 mIU/mL in urine. The test utilizes a combination of monoclonal and polyclonal antibodies to selectively detect elevated levels of hCG in urine. At the level of claimed sensitivity, FRENOVO Human Chorionic Gonadotropin (hCG) Rapid Test shows no cross-reactivity interference from the structurally related glycoprotein hormones hFSH, hLH and hTSH at high physiological levels.

The test utilizes a combination of antibodies including mouse monoclonal anti-hCG antibodies to selectively detect elevated levels of hCG. The assay is conducted by adding urine specimen to the specimen well of the test cassette and observing the formation of colored lines. The specimen migrates via capillary action along the membrane to react with the colored conjugate. Positive specimens react with the specific colored antibody conjugates and form a colored line at the test line region of the membrane. Absence of this colored line suggests a negative result. To serve as a procedural control, a colored line will always appear at the control line region if the test has been performed properly.

### MATERIALS PROVIDED

Each kit contains:

1. Test Devices: 20 pieces test devices individually pouched.
2. Droppers: 20 pieces droppers of 25 ul
3. Package insert: 1 piece attached.

### MATERIALS REQUIRED BUT NOT PROVIDED

- Timer or stopwatch.
- Specimen collection containers
- Disposable gloves and/or protective clothing
- Micropipette

### WARNINGS

1. Read the package insert completely before using the product. The instructions must be followed carefully as not doing so may result in inaccurate results.
2. FRENOVO Human Chorionic Gonadotropin (hCG) Rapid Test is for diagnostic use only.

FRENOVO Human Chorionic Gonadotropin (hCG) Rapid Test is for professional use only.

3. Perform test at room temperature.

### PRECAUTIONS

1. FRENOVO Human Chorionic Gonadotropin (hCG) Rapid Test is for professional use only.
2. The package insert instructions must be followed to ensure optimum test performance.
3. The used test device should be discarded according to local regulations.

### Handling Precautions

1. Do not use if the kit box safety seal is absent, damaged or broken.
2. Do not use any device if the pouches have been perforated.
3. Each device is for single use only.
4. Do not use the kit past the expiration date (this date is printed on the kit box).
5. Adequate lighting is required to read the test results.
6. The result should be read immediately after the end of the 5 minutes incubation time following the addition of specimen and wash buffer solution. Do not read results beyond 10 minutes.

### STORAGE INSTRUCTIONS

1. FRENOVO Human Chorionic Gonadotropin (hCG) Rapid Testkit should be stored between 2-30°C and the shelf life is 36 months.
2. If stored refrigerated, ensure that the pouched device is brought to room temperature before opening.
3. Do not freeze the kit.

### SAMPLE COLLECTION AND PREPARATION

The urine specimen must be collected in a clean and dry container. Urine collected at any time of the day may be used. Urine specimens exhibiting visible precipitates should be centrifuged, filtered, or allowed to settle to obtain a clear supernatant for testing. Urine specimens may be stored at 2-8° C for up to 48 hours prior to testing. For prolonged storage, specimens may be frozen and stored below -20° C. Frozen specimens should be thawed and mixed well before testing.


### QUALITY CONTROL

An internal procedural control is included in the test. a colored line appearing in the control line region (C) is an internal valid procedural control, it confirming adequate membrane wicking. Control standards are not supplied with this kit; however, it is recommended that positive and negative controls be tested as a good laboratory practice to confirm the test procedure and to verify proper test performance.

### TEST PROCEDURE

**Allow the test, urine specimen, and/or controls to reach room temperature (15-30°C) prior to testing.**

1. Bring the pouch to room temperature before opening it. Remove the test device from the sealed pouch and use it as soon as possible.
2. Place the test device on a clean and level surface. Hold the dropper vertically and transfer 2 full drops of urine (approx. 50 µL) to the specimen well of the test device, and then start the timer. Avoid trapping air bubbles in the specimen well. See the illustration below.
3. Wait for the colored line(s) to appear. Read results at 5 minutes. Do not interpret the result after 10 minutes.



### INTERPRETATION OF RESULTS

(Please refer to the illustration above)

**Negative result:** if there is only a quality control line C, the detection line is colorless, indicating that hCG has not been detected and the result is negative.

**Positive result:** if both the quality control line C and the detection line T appear, HCG has been detected and the result is positive.

**Invalid result:** if the quality control line C is not observed, it will be invalid regardless of whether there is detection line, and the test shall be conducted again.

#### LIMITATIONS

1. Very dilute urine specimens, as indicated by a low specific gravity, may not contain representative levels of hCG. If pregnancy is still suspected, a first morning urine specimen should be collected 48 hours later and tested.
2. False negative results may occur when the levels of hCG are below the sensitivity level of the test. When pregnancy is still suspected, a first morning urine specimen should be collected 48 hours later and tested.
3. Very low levels of hCG (less than 50 mIU/mL) are present in urine specimen shortly after implantation. However, because a significant number of first trimester pregnancies terminate for natural reasons, a test result that is weakly positive should be confirmed by retesting with a first morning urine specimen collected 48 hours later.
4. A number of conditions other than pregnancy, including trophoblastic disease and certain non-trophoblastic neoplasms including testicular tumors, prostate cancer, breast cancer, and lung cancer, cause elevated levels of hCG. Therefore, the presence of hCG in urine specimen should not be used to diagnose pregnancy unless these conditions have been ruled out.
5. This test provides a presumptive diagnosis for pregnancy. A confirmed pregnancy diagnosis should only be made by a physician after all clinical and laboratory findings have been evaluated.

#### PERFORMANCE CHARACTERISTICS

##### Accuracy

A clinical evaluation was conducted comparing the results obtained using the FRENOVO Human Chorionic Gonadotropin (hCG) Rapid Test and another commercially available urine hCG test. The study included 200 urine specimens and both assays identified 100 negative and 100 positive results. The results demonstrated 100% overall agreement (for an accuracy of > 99%) of the FRENOVO Human Chorionic Gonadotropin (hCG) Rapid Test.

##### Cross-Reactivity

A study was conducted to determine the cross-reactivity of the test with compounds in either negative urine or positive urine. The test has been standardized to the W.H.O. Fourth International Standard (75/589). The addition of LH (300 mIU/mL), FSH (1,000 mIU/mL), and TSH (1,000  $\mu$ IU/mL) to negative (0 mIU/mL hCG) and positive (20 mIU/mL hCG) specimens showed no cross-reactivity with FRENOVO Human Chorionic Gonadotropin (hCG) Rapid Test.















##### Interfering substances:

A study was conducted to determine the following potentially interfering substances, which were added to hCG negative and positive specimens, then tested with FRENOVO Human Chorionic Gonadotropin (hCG) Rapid Test, tested no interfered.

Compounds	Conc.mg/dL	Compounds	Conc.mg/dL
Acetaminophen	20	Acetone	1000
Acetylsalicylic Acid	20	Acetoacetic Acid	2,000
Ampicillin	20	Ascorbic Acid	20
Atropine	20	Albumin	2,000
$\beta$ -Hydroxybutyrate salt	2000	Benzoylcegonine	10
Bilirubin	20	Brompheniramine	20
Caffeine	20	Cannabinol	10
Clomiphene	100	Cocaine	10
Codeine	10	Cholesterol	500
Creatine	20	Dextromethorphan	20
DMSO	5%	EDTA	80
Ephedrine	20	Ethanol	1%
Estril	2	Estrone 3-Sulfate	10
Gentisic Acid	20	Glucose	2,000

Hemoglobin	1,000	Heroin	1
Ibuprofen	20	Methadone	10
Methamphetamine	10	Methanol	10%
Morphine	0.6	Oxalic Acid	40
Phenothiazine	20	Phenylpropanolamine	20
Pregnanediol	2	Salicylic Acid	20
Tetracycline	20	Triglycerides	1,200
Theophylline	20	Urea	2,000
Uric Acid	20	-----	

#### INDEX OF SYBOML

	In vitro diagnostic medical device		single-use, Please don't reuse it
	Use-by date		Consult instructions for use
	Cautions		Manufacturer
	Temperature limit		Batch code
	Date of manufacture		Keep Dry
	Avoid overexposure to the sun		Don't use the product when the package is damaged
	CE mark		Biological risks



**Hangzhou Frenovo Biotech Co., Ltd.**  
**Address:** Room 401, Building 36, No. 488-1, Donghu North Road, Donghu Street, Yuhang District, Hangzhou City, Zhejiang Province, China.  
**Tel:** 86-0571-89170657  
**Email:** business@frenovo.com



**Lotus NL B.V.**  
**Address:** Koningin Julianaplein 10, 1e Verd, 2595AA, The Hague, Netherlands.  
**Email:** peter@lotusnl.com

#### INSTRUCTION APPROVAL AND REVISION DATE

Approval Date:  
Revision Date:  
Date of Issue: