



## Drug Of Abuse Test ▼

Product Name	Specimen	Format	Certificate
Morphine (MOP) Rapid Test	Urine	Cassette/Sheet	CE
Cocaine (COC) Rapid Test	Urine	Cassette/Sheet	CE
Amphetamine (AMP) Rapid Test Kit	Urine	Cassette/Sheet	CE
Barbiturate (BAR) Rapid Test Kit	Urine	Cassette/Sheet	CE
Benzodiazepines (BZO) Rapid Test Kit	Urine	Cassette/Sheet	CE
Ecstasy (MDMA) Rapid Test Kit	Urine	Cassette/Sheet	CE
Marijuana (THC) Rapid Test Kit	Urine	Cassette/Sheet	CE
Methadone (MTD) Rapid Test Kit	Urine	Cassette/Sheet	CE
Methamphetamine (MET) Rapid Test Kit	Urine	Cassette/Sheet	CE
Mutli-Drugs (2-10 Strips) Rapid Test Cup	Urine	Cassette/Sheet	CE
Mutli-Drugs (2-10 Strips) Rapid Test Panel	Urine	Cassette/Sheet	CE
Opiate (OPI) Rapid Test Kit	Urine	Cassette/Sheet	CE
Tricyclic Antidepressants (TCA) Rapid Test Kit	Urine	Cassette/Sheet	CE

## Immunofluorescence-Platform ▼

Product Name	Working range	Performance	LOD	Assay time	Sample vol.
cTnI	0.1-40ng/mL	CV≤15%	0.1ng/mL	15mins	100ul
CRP	0.5mg/L~200mg/L	CV≤15%	0.5mg/L	3mins	100ul
PCT	0.1-100ng/mL	CV≤15%	0.1ng/L	15mins	100ul
β-HCG	5~20000mIU/mL	CV≤15%	5mIU/mL	15mins	100ul
HbA1c	3%~14%	CV≤15%	3%	5mins	100ul



- Website -



- Catalogue -